Bioliss is a game-changing **premium product** designed with professional hairdressers in mind. Its advanced formula, enriched with the natural power of **organic acids**, not only delivers flawlessly straightened hair but also ensures results that **last for up to 8 months**.*

What sets Bioliss apart is its **unique recipe** that preserves the hair's disulphide bonds while strengthening its structure, leaving hair more resilient to everyday challenges. And the best part? Bioliss contains **no harmful ingredients**, making it a safe and gentle choice for all hair types. Bioliss combines cutting-edge innovation with care, giving your clients stunning results they'll love.

BEAUTY & CREATIVI



What's in the Bioliss treatment kit? Clarifying Shampoo 100ml or 500ml and Nourishing and Straightening Fluid 250ml or 1000ml.

Bioliss Bioplasty is more than just a straightening treatment – it also deeply hydrates and improves hair porosity. With Bioliss, you'll achieve beautifully straight hair on any curl pattern, while creating soft, silky, and glossy strands. It also guarantees the complete elimination of frizz, leaving hair looking flawless.

Bioliss Bioplasty is a hybrid treatment that's not only highly effective and long-lasting but also incredibly versatile. Thanks to its **perfectly balanced pigments**, this treatment works beautifully **on all hair types** – and is especially suited for blonde hair!

Whether your client's hair is natural, bleached, or coloured, Bioliss delivers **exceptional results without compromising the health of the hair**.

Bioliss Bioplasty works like a **non-invasive hair lightening treatment**. Coloured hair can lighten by up to 3 levels, revealing a stunning, cool blonde tone, while natural hair remains unaffected, preserving its original shade.

This makes Bioliss an ideal choice for enhancing and refining blonde tones, all while maintaining the integrity of the hair.



Bioliss Bioplasty is free from formaldehyde, carbocysteine, and their derivatives. It's also 100% vegan, offering a safe and ethical choice for salons and clients alike. The Bioliss Bioplasty treatment is **thermoactive**,

meaning that for the best results, the hair should be blow-dried after each wash while being gently brushed or combed. This step helps to fully activate the treatment and maintain its stunning effects.

A simple routine for your clients to keep their hair looking its absolute best!

Step 1: Bioliss Shampoo

1. Wash the hair with the Bioliss clarifying shampoo included in the kit to remove any product build-up, such as serums or conditioners.

After rinsing, towel-dry the hair to remove excess water and dry it, leaving the appropriate level of moisture suited to the hair type and condition.

Moisture levels before product application:

0% - Virgin or Afro hair: Apply the product to completely dry hair.

10% - Hair in good condition: Leave 10% moisture. If the hair begins to dry during the treatment, rehydrate the strands with water using a spray bottle.

20% - Bleached hair: Leave 20% moisture. Rehydrate as needed during the process with a spray bottle.

Step 2: Bioliss Fluid Application

1. Apply the Product: Apply Bioliss Fluid 1cm away from the roots, evenly covering the entire length of the hair – avoid applying directly to the scalp. For Afro-textured hair, apply close to the roots, but first perform a patch test to check the scalp's reaction. If there's any sensitivity, maintain a 1–2cm gap from the roots. Ensure the product is applied generously and evenly, covering each strand thoroughly. It's better to use slightly more than too little.

2. No Rinsing Required: Bioliss Fluid does not require rinsing. After the processing time, gently blot the hair with a towel to remove excess product. In this case, the hair must be washed in the salon immediately after the treatment – no need to wait 48 hours for the first wash. Exception: For Afro-textured hair, wait 48 hours before the first wash.

3. Completely Dry the Hair: This step is crucial. Never straighten damp or wet hair, as it can cause permanent damage.

Divide the hair into four sections and dry thoroughly, strand by strand, from roots to ends using a brush. This stage activates the product fully and ensures longer-lasting, beautiful results.

4. Straightening: Begin by dividing the hair into four sections. Straighten very thin strands using a titanium flat iron with wide plates. Slowly and carefully press each strand until a glossy finish is achieved.

Adjust the temperature and number of passes according to the hair type and condition, following the table below.

hair type	flat iron temperature (°C)	number of passes	processing time (min)
bleached, delicate, fine, high-porosity, straight/frizzy	170	14 - 18	40
bleached and thick, wavy or slightly curly	180 - 200	18 - 24	50
medium-porosity, moderately sensitised, frizzy, wavy/curly, coloured, or virgin	200-215 coloured hair 210-220 virgin hair	18-20 coloured hair 22-27 virgin hair	55
thick, strong curls, coily, virgin, -textured-Afro	220-225	24-30	70 - 90

Bioliss Bioplasty Treatment

Additional Recommendations

• Before proceeding, perform a **wet strand test** to determine if the hair is suitable for straightening. **If the hair stretches, frays, or breaks under tension, the treatment should be avoided.** In such cases, recommend trimming the damaged areas and following up with intensive hydrating and restorative treatments, such as **Platinum B-tox**.

• For clients with sensitive skin or allergy concerns, **conduct a patch test 24 hours before the treatment**. A mild tingling sensation during the procedure is normal due to the presence of acids in the formula.

• During the treatment, ensure that **metal tools such as combs, clips, or bowls are not used**, as they may interfere with the product's performance.

Contraindications: Active scalp conditions, allergies to the product's ingredients, or severely damaged hair.

For the Bioliss treatment, it is essential to use a titanium straightener with wide plates (approximately 4.5 cm). The straightener must maintain a consistent temperature throughout the procedure. When choosing a straightener, ensure it allows temperature adjustments in increments of at least 5°C for optimal precision.

. Our Diana Titanium Straightener is the perfect tool for this treatment and is available for purchase on our website.

Before proceeding, ensure the client fully understands the treatment process. **Provide them with a leaflet**, and encourage them to **sign a client form** (available on our website). This form should note key details, such as the product mix and the number of passes with the straightener, as these are valuable for follow-up services.

For **sensitive or dry hair**, we recommend preparing client's hair with products like **Pre Nano B-tox**, which hydrates and softens, and **Nanoplex**, which protects the hair structure from the straightener's heat.



First Hair Wash After Treatment

Hair should be washed immediately after the treatment (allow 5–10 minutes for the hair to cool down after straightening). We recommend using the dedicated **HomeCare Bioliss shampoo and conditioner**. Alternatively, you can use the Nano shampoo and mask from the Nanoplasty line.

For thick, curly, or Afro-textured hair, it is advised to wait 48 hours before the first wash.

Steps for washing:

1. Wash the hair thoroughly with **Bioliss shampoo**, ideally twice, to ensure the product is fully removed from the hair.

2. Towel-dry the hair to remove excess water.

3. Apply Bioliss conditioner evenly throughout the hair.

4. Leave the conditioner on for 5–10 minutes. During this time, use a wide-tooth comb to gently detangle the hair.

5. Rinse the conditioner out thoroughly with warm or lukewarm water (avoid using very hot water).

Post-Treatment Care with Bioliss

The Bioliss range includes aftercare products specifically designed to maintain and extend the results of the treatment.

BIOLISS SHAMPOO

DIANY

BIOLISS

MAINTENANCE

HOMECARE TREATMENT

250 ml

This gentle formula is kind to both the scalp and hair strands, preserving the scalp's hydrolipid barrier and protecting the hair along its length. The shampoo optimises scalp hydration while leaving hair moisturised, soft, and shiny.

BIOLISS MASK

The mask combines hydrating, emollient, and protein-rich ingredients, offering comprehensive protection after straightening treatments. It helps to prolong the effects of the treatment, keeping hair healthy and vibrant.



The formulas of Bioliss aftercare products are specifically designed to complement the ingredients used in the straightening treatments. They also feature an optimised pH level to help clients enjoy the results of their treatment for as long as possible.

Using shampoos or conditioners with harsh or incompatible ingredients, including clarifying or drugstore products, can significantly reduce the longevity of the treatment's effects.

Order Bioliss products at www.dianacosmetics.nl

